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RESORBABLE COLLAGEN MEMBRANE

ENG INSTRUCTIONS FOR USE

Resorbable collagen membrane is an implantable material that is indicated for guided tissue regeneration (GTR) procedure in periodontal defects.

DESCRIPTION

Resorbable collagen membrane is from purified collagen fibres type I and ΙΙΙ derived from skin porcine origin.

Resorbable collagen membrane is fully resorbable eliminating the need for a second surgical procedure often required to remove a non-resorbable membrane.

Resorbable collagen membrane protects the implantation site from invasion of soft tissues (epithelial and conjunctival cells) and ensures this barrier function during 12 weeks. This membrane allows the passage of biological fluids and is intended to be resorbed in 24 weeks (lifetime).

The biocompatibility of Resorbable collagen membrane has been demonstrated.

INDICATIONS

Resorbable collagen membrane is intended for use in periodontal / dental surgery procedures as a barrier to epithelio-conjunctive colonization in the area of periodontal defect, dental implant, bone defect or ridge reconstruction.

Considering indications, the combination with a bone graft may help bone growth for guided bone regeneration (GBR).

CONTRAINDICATIONS

- The Resorbable collagen membrane should not be placed on the site of an active infection.

- Refusal of the patient to adhere to surgical follow-up and to the limit in activity level.

- Fever and/or local inflammation.

- It should only be used in the surgeries indicated.

PRECAUTIONS FOR USE

The content of the double pouches is designed for single use only.

Any salivary or other contamination of the material and surgical site should be minimized in order to avoid bacterial contamination.

Before laying, the surgeon should ensure that active or recent infections have been treated properly.

If endosseous implants are involved, the membrane should be used only in combination with a stable implant and not in lieu of achieving primary implant stability.

The smallest movement on the tissue underneath need to be avoided.

The material has not been tested on pregnant women or for pediatric use.

RECOMMENDED IMPLANTATION TECHNIQUE

It is advisable that only clinicians trained in related treatment planning and in the technique of placing periodontal membranes should use Resorbable collagen membrane.

Resorbable collagen membrane should be removed from the packaging using sterile gloves or instruments.

Mucoperiosteal incision flaps are developed in the site to be treated. The incision should be sulcular as possible. Clinician should perform thorough debridment. As much as possible tissue must be preserved to allow for primary closure of the wound and correct positioning of the flaps.

If necessary, bone defect could be filled with bone grafting material.

Resorbable collagen membrane can be cut with sterile scissors to desire shape.

Resorbable collagen membrane can be placed either dry or hydrated. For better handling characteristics, the membrane can be hydrated in sterile water or saline solution prior to the final placement especially when the practitioner intends to conform or bend the membrane at the site of implantation.

Resorbable collagen membrane should completely cover the defect and provide a minimum 2 mm bone-to-membrane contact along the periphery of the defect to prevent invasion by soft tissues.

The membrane can be fixed by suturing. This can be indicated to prevent the membrane from moving.

The mucoperiosteal flap is sutured over the collagen membrane (e.g., single sutures and deep mattress sutures).

The wound should be closed completely to avoid exposure of the membrane to the open air and the oral environment.

Special instructions for use in periodontology

A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planning, and occlusal adjustment when indicated.

A postoperative maintenance phase can help to ensure long-term therapeutic success.

In order to avoid the formation of a long junctional epithelium, Resorbable collagen membrane must be adapted closely to the tooth (e.g. with additional fixation using suture material).

Post-operative care

If the membrane becomes exposed, the dehiscence usually heals by itself within several weeks. Membrane removal is usually not necessary. However, to minimize bacterial contamination, rinsing with bactericidal solutions is recommended.

In the event that a membrane removal is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curetted to remove any inflamed or infected tissue.

ADVERSE EFFECTS

Since Resorbable collagen membrane is a collagen product allergic reactions may not be totally excluded and adverse reactions like fever experienced with another microfibrillar collagen may be related.

The barrier function of the membrane is no longer assured if the resorption of the membrane is faster than the regeneration of the tissues.

STERILISATION AND STORAGE

Do not use after expiry date.

Resorbable collagen membrane is Bêta irradiated sterilized in double aluminium pouches.

Cannot be re-sterilized. The device is sterile unless the package has been opened, damaged, or otherwise contaminated.

The resorbable collagen membrane is sold sterile and can not be reused once the package is opened. The consequences of a re-sterilisation and/or a re-use of the Resorbable collagen membrane have not been evaluated during the design process. So, the mechanical and resorption performances of Resorbable collagen membrane are not guaranteed in case of re-sterilisation and/or re-use.

This product must be handled by suitably trained and qualified personnel having read and understood this instruction for use. Store at room temperature 15-25 °C (59-77 °F) and in a dry place.

SYMBLES	Latest revision date: 07/2017	
 Single use	 Expiration date	 Quantity
 This product must not be re-sterilized	 Consult instructions for use	 Temperature limit 15°C/25°C (590-770°F)
 Do not use if package is damaged	 Sterilization by irradiation	
 Manufacture date (Sterilization)	 Manufacturer	

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