**+ADM-1st Choice Medicals**  
smart & quality solutions at a reasonable price

# Ora-Aid

A New Conceptual Intraoral Wound Dressing Strip

## Ora-Aid Application


Extraction

Socket grafting

Free gingival graft

Temporary Crown

"Ora-Aid Video"

**TBM**

TBM Corporation

ISO  
13485

# Extraction (OAC Case)

**Reference information >>**  
Please refer to our QR code for various clinical information and testimonials

- Acute Oroantral communication (OAC), a severe case caused by the extraction of a maxillary (upper) tooth (usually a maxillary first molar), has usually been treated with SG or PRF due to the lack of a reasonable and efficient treatment method. Ora-Aid is an innovative solution for this special case of extraction as well as for general extraction cases, allowing for the closure of OAC in a less invasive and more practical way.

Article: Is an Attachable Oral Wound Dressing effective at closing an Acute Oroantral Communication(OAC)?  
- Yusuf Nuri Kaba, DDS,\* Emrah Soylu, DDS, PhD, Ahmet Emin Demirbas, DDS, PhD and Musab Suleyman Kilavuz, DDS

| In OAC Case          | Ora-Aid  | Sterile Gauze   | PRF   |
|----------------------|--|---|---|
| <b>Advantages</b>    | Reasonable price, same recovery outcomes compared to other methods, but with the lowest VAS scores and dramatically shorter surgery time due to easy application and easy-to-handle materials. | The most common techniques for protecting against blood clots that may form in the socket.  | PRF accelerates soft tissue recovery by increasing fibroblast proliferation.  |
| <b>Disadvantages</b> | None   | The porous gauze structure increases the risk of secondary infections, plaque buildup, and unpleasant odors. It's also a <b>hassle to remove</b> the gauze afterward. | Requiring additional technical equipment, preparation time. Average 10~15min). Additionally, collecting blood from patients may cause anxiety, and the preparation process may increase contamination risk. |

|                         |   |  |  |
|-------------------------|---|--|--|
| <b>Ora-Aid Features</b> | <ul style="list-style-type: none"> <li>Ora-Aid helps <b>hemostasis</b> of the surgical site and rapid <b>natural healing</b> of soft tissues, dramatically <b>reducing patient pain and discomfort</b></li> <li>Flexible and self-adhesive for <b>easy manipulation</b>.</li> <li>Quickly adheres to surfaces for hassle-free placement.</li> </ul> |  |  |
| <b>Clinical Case</b>    |  <p>&lt; Application of Ora-Aid after extraction of palatal lateral incisors &gt;</p>   |  |  |

# Socket grafting

**Reference information >>**  
Please refer to our QR code for various clinical information and testimonials

Acute Oroantral < communication(OAC) closing with Ora-Aid>

| Product                | Membrane   | Ora-Aid  | Remark   |
|------------------------|--|--|--|
| <b>Price</b>           | Expensive  | Very Reasonable  | Some membrane price is more than US\$100.00 for 1ea  |
| <b>Sequence</b>        | Place Membrane inside of flap: Not easy<br>Close flap: May need released incision and apply tension to closeflap<br><b>and it may damage to attached ligament. It is a very difficult procedure.</b> | If socket size is less than 12mm:<br>Attach Ora-Aid on the wound and suture it according to provided protocol<br>If socket size is bigger than 12mm:<br>Stitch on the wound and same procedure with above. | <b>When wall defected, it should be use Membrane.</b><br>Ora-Aid should always be placed over the wound and sutured according to our protocol. |
| <b>Clinical Case</b>   |  <p>&lt;Bond Apatite + Ora-Aid + Suture&gt;</p>  |  |  |
| <b>Clinical Result</b> | No difference  |  | Ora-Aid is more economical and easier to use.  |
| <b>Usability</b>       | Absorbable Membranes:<br>Relatively good maneuverability, but less space for bone regeneration.<br>Non-absorbable Membranes:<br>Relatively difficult to handle and place.                            | Easy to trim in accordance with wound size, shape and to apply.<br>Quickly adheres to the surface, saving time with no hassle.   | A simple figure 8 suture can be placed over the Ora-Aid to provide longer protection for the area.   |

## FGG (Free Gingival Graft)

Reference information >>

Please refer to our QR code for various clinical information and testimonials

How FGG cases have been treated so far:

- Suturing: Just painful.
- Collagen tape and sutures: simple, but costly, with little pain relief.
- 3D printing: Requires scanning, 3D printing, and modification; The process is cumbersome and expensive, patients tend to be uncomfortable wearing it.
- Flowable resin composite stent: uncomfortable, not effective in reducing pain.

Ora-Aid : An innovative solution for FGG cases that significantly reduces pain associated with the donor site, eliminates discomfort and foreign bodysensation, and creates the best environment for healing, enabling faster hemostasis and healing.

### Ora-Aid

|               |   | Remark   |
|---------------|---|--|
| Price         | More than reasonable  | Ora-Aid: 20ea per box  |
| Sequence      | After suturing the donor site, simply attach Ora-Aid above the suture.<br>We recommend suture if it needs to protect wounds around week.  | Suture should be done by provided suture protocol.   |
| Clinical Case |    <br><Ora-Aid + Suture>                      <Surgical + Ora-Aid + Suture> |  |
| Usability     | The most effective way to <b>relieve pain</b> and support <b>healing procedure</b> yet easy to apply.<br>Quickly adheres to surfaces for hassle-free.   | Collagen tape at the donor site or using a flowable resin after applying Ora-Aid could be even <b>more effective</b> . |

## Temporary Crown

Reference information >>

Please refer to our QR code for various clinical information and testimonials

Ora-Aid is an innovative solution that is significantly effective against hemostasis. It reduces pain from external stimuli after a temporary crown, eliminates discomfort and foreign body sensation that patients may feel from mild gum pain and discomfort: After a temporary crown is placed, it is common to feel mild pain or discomfort around the gum line due to the artificial manipulation of the gum tissue. This may manifest as a dull ache, sensitivity to hot or cold foods, or slight throbbing. Also, Ora-Aid creates the ideal environment for healing and support through the entire healing process.

### Ora-Aid

|               | Detail   | Remark   |
|---------------|--|--|
| Price         | More than reasonable   | Ora-Aid: 20ea per box  |
| Sequence      | After suturing the surgical site, place an Ora-Aid, and attach a temporary crown over it.  | Temporary crown anchors the Ora-Aid and holds it for an extended period of time(case1) without requiring a suture. |
| Clinical Case | <div>■ case 1</div>   <div>Soft tissue after 4weeks, Ora-Aid is still attached</div> <div>■ case 2</div>    <div>Case with Ora-Aid attached before temporary denture</div> |  |
| Features      | <ul style="list-style-type: none"><li>■ Ora-Aid helps <b>hemostasis</b> of the surgical site and rapid <b>natural healing</b> of soft tissues, dramatically <b>reducing</b> patient <b>pain and discomfort</b>.</li><li>■ Flexible and self-adhesive for <b>easy manipulation</b>.</li><li>■ Quickly adheres to surfaces for hassle-free.</li></ul>  |  |

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+ADM-1st Choice Medicals  
 smart & quality solutions at a reasonable price

Web: [www.fcm.com.gr](http://www.fcm.com.gr)  
 Email: [info@fcm.com.gr](mailto:info@fcm.com.gr)  
 Tel.: +30 211 42 13 059  
 Cyprus: +357 97740827