



www.lazonlaser.com

<https://youtu.be/teuKagI4JpE> (promo video)

ADVANTAGES:

- > 1. OPERATION ON IPHONE / ANDROID, WHICH MAKES IT A GREAT EXPERIENCE.
- > 2. BEST CONFIGURATIONS IN THE MARKET.
- > 3. EXCELLENT SUPPLIERS FROM GERMANY, DENMARK, US AND ENGLAND.
- > 4. COMPETITIVE PRICE IN THE MARKET.
- > 5. 2-YEAR WARRANTY PERIOD.
- > 6. CE APPROVED.
- > 7. WIRELESS
- > 8. CHOOSE BETWEEN TWO WAVELENGTHS 810 Nm & 980 Nm
- > 9. THE STRONGEST IN POWER, IN COMPARISON WITH OTHER DIODE LASER (see comparative chart)
- > 10. BETTER POWER STABILITY (see comparative chart)
- > 11. BETTER ACCURACY chart
- > 12. BETTER COOLING METHOD (see comparative chart)
- *** (the best in 25 categories – see comparative chart)
- > 13. CHOOSE BETWEEN DIFFERENT TYPES OF HANDPIECE
- > 14. 3 TYPES OF TIPS (200 µm – 300 µm – 400 µm) + bending tool
- > 15. THE MOST POWERFUL BATTERY
- > 16. VERY SMALL & LIGHT
- > 17. VERY SMART IOS/ ANDROID APPLICATION (which even has patient records inside with pictures or videos)

-
- laser for soft tissue operation also teeth whitening and
 - bio stimulation with competitive price and unique design. In the near
 - future, laser treatment will be soon suitable for almost everyday practice.
 - And Solase has got widely recognition from distributors and dentists
 - on three exhibitions from March this year.

MARKETING VIDEOS:

https://youtu.be/3tNP_iFi7nQ (quick start guide)

<https://youtu.be/27m7CPK5YRY> (quick start guide)

<https://youtu.be/d9TiwoNXHyU> (instructions)

<https://www.lazonlaser.com/product/jiguangqi.aspx>

<https://youtu.be/teuKagI4JpE>

<https://youtu.be/KZYuzsIXXew> (accuracy and stability test)

<https://youtu.be/0ENwiw9BbUw> (precautions)

<https://youtu.be/d9TiwoNXHyU> (application introduction)

CLINICAL VIDEOS:

<https://youtu.be/jlxH0DKxLtk> (herpes)

<https://youtu.be/FdlSg-eoESU> (gingiva pigmentation)

<https://youtu.be/KuRs4dP3plM> (periodontal pocket)

<https://youtu.be/GxkKRQJjwdM> (periodontal pocket disinfection)

TECHNICAL SPECIFICATIONS:

Wavelength - Two models, 808nm or 976nm

Power Output -

-976nm: 0.1-16W@pulse mode, 0.1-10W@CW mode

-808nm: 0.1-12W@pulse mode, 0.1-7W@CW mode

Pulse Duration - 10μs-0.9s

Pulse Repetition Rate - Up to 50kHz

Tip Diameter - 200μm, 300μm, 400μm

Power Supply Input -100-240VAC, 50-60Hz

Battery Pack - Li-ion Rechargeable, Last for one week with normal use

Size / Weight - 138 x 124 x 168 mm, 1.3 kg

APPLICATIONS:

Surgery

Gingivectomy/Gingivoplasty

Implant Uncovery

Operculectomy

Frenectomy/Frenotomy

Hemostasis

Troughing

Aphthous Ulcers

Fibroma

Epulis

Crown Lengthening

Periodontology

Periodontal Pockets

Peri-implantitis

Endodontics

Pulpotomy

Root Canal Disinfection

Bio-Stimulation

Temporary Pain Relief

TMJ Disorders

Dentine Hypersensitivity

Cosmetics

Tooth Whitening

Gingival Recontouring

Gingival Pigmentation