

BLDC MOTOR SERIES

OZ

60K 8.0 N·cm speed

- Digital RPM Display
- Digital watch
- Strong torque at low speed
- Vertical / horizontal placement
- Hand ↔ Foot control select
- Overload protection system
- Easy carry with handle
- Slim & Lightweight handpiece



OZ100

OZ Elite

50K 8.2 N·cm speed

- Digital RPM Display
- Strong torque at low speed
- Vertical / Horizontal placement
- Hand ↔ Foot control select
- Overload protection system
- Easy carry with handle
- Slim & Lightweight handpiece



H180

TRAUS B110

50K 8.2 N·cm speed knee

- Digital RPM Display
- Knee ↔ Foot control select
- Compact design
- Auto-Cruise function
- Overload protection system
- Slim & Lightweight handpiece



H180

OZ Black

50K 7.8 N·cm speed

- Digital RPM Display
- Excellent durability
- Vertical / Horizontal placement
- Hand ↔ Foot control select
- Overload protection system
- Easy carry with handle



H500

Forte 200a

50K 8.2 N·cm speed knee

- Digital RPM Display
- Knee ↔ Foot control select
- Ergonomically designed
- Excellent durability
- Auto-Cruise function



H500

Forte 400P

50K 7.8 N·cm speed

- Digital RPM Display
- Desktop type
- Hand ↔ Foot control select
- Excellent durability
- Overload protection system



FORTE 100IIP

TRAUS B150

50K 7.8 N·cm speed

- Digital RPM Display
- Desktop type
- Compact-Stylish Body
- Hand ↔ Foot control select
- Excellent durability
- Overload protection system



FORTE 100IIP



OZ100

H180

H500

F100IIP

F100aEI



SAESHIN

T. 82-53-587-2341 | F. 82-53-580-0999 | E. info@saeshin.com

www.saeshin.com



DC MOTOR SERIES

STRONG B140D

45K 4.0 N·cm speed

- Digital RPM Display
- Powerful torque
- Hand ↔ Foot control select
- Excellent durability
- Overload protection system
- 10 drill bits holder on control box
- Variable speed control pedal



H450

STRONG 210

35K 3.0 N·cm on/off speed

- Built-in handpiece holder
- 16 drill bits holder on control box
- Hand ↔ Foot control select
- Overload protection system
- Variable speed control pedal (optional)



STRONG 105L

Brilliant B350D

35K 3.0 N·cm on/off

- Digital RPM Display
- Built-in handpiece holder
- 8 drill bits holder on control box
- Hand ↔ Foot control select
- Overload protection system



H350

Brilliant B350

35K 3.0 N·cm on/off

- Built-in handpiece holder
- 8 drill bits holder on control box
- Hand ↔ Foot control select
- Overload protection system



H350

Brilliant 210

35K 3.0 N·cm on/off

- Built-in handpiece holder
- 16 drill bits holder on control box
- Hand ↔ Foot control select
- Overload protection system



H350



H450

H400

H350

105

108E

PORTABLE SERIES

Brilliant B170

30K 3.0 N·cm on/off

- Powerful torque with BLDC motor
- Digital RPM Display
- Slim & Lightweight handpiece
- Battery 10 hours use, with 2 hours charging time
- One-click pause function
- Hand ↔ Foot control select



H175

H170E

THUMB-D

30K 3.0 N·cm

- Digital RPM Display
- Slim & Lightweight handpiece
- Battery 8 hours use, with 1 hours 30 minutes charging time
- One-click pause function



H325

PEHP